



CIE AS Chemistry



Your notes

22.1 Infrared Spectroscopy

Contents

* Infrared Spectroscopy



Your notes

Infrared Spectroscopy

Interpreting IR Spectra

- **Infrared (IR) spectroscopy** is a technique used to identify compounds based on changes in **vibrations** of atoms when they **absorb** IR of certain **frequencies**
- A **spectrophotometer** irradiates the sample with electromagnetic waves in the infrared region and then detects the **intensity** of the wavelength of IR radiation which goes through the sample
- All organic molecules absorb IR radiation and depending on which energies of radiation are absorbed, bonds between atoms will vibrate by **stretching, bending and twisting**
- The molecules will only vibrate at a specific frequency
- The **resonance frequency** is the specific frequency at which the molecules will vibrate to stimulate larger vibrations
- Depending on the rest of the molecule, each vibration will absorb specific wavelengths of IR radiation which are also shown as the **reciprocal** of the wavelength
 - This unit is called the **wavenumber** (cm^{-1})
- Particular absorbances have characteristic **widths** (broad or sharp) and **intensities** (strong or weak)
 - For example, hydrogen bonds cause the O-H bonds in alcohols and carboxylic acids to be **broad** whereas the C=O bond in carbonyl ($\text{C}=\text{O}$) groups have a strong, sharp absorbance peak
- The energies absorbed by different functional groups are given as a range and an unknown compound can be identified by comparing its IR spectrum to the IR spectrum of a known compound

Absorption range of bonds table

Bond	Functional groups containing the bond	Characteristic infrared absorption range (in wavenumbers) / cm^{-1}
C-O	Hydroxy, ester	1040 - 1300
C-C	Aromatic compound, alkene	1500 - 1680
C=O	Amide	1640 - 1690
	Carbonyl, carboxyl	1670 - 1740
	Ester	1710 - 1750
$\text{C} \equiv \text{N}$	Nitrile	2200 - 2250
C-H	Alkane	2850 - 2950
N-H	Amine, amide	3300 - 3500

O-H	Carboxyl	2500 - 3000
	Hydroxyl	3200 - 3600

- Due to some absorption bands overlapping each other, other analytical techniques such as mass spectroscopy should be used alongside IR spectroscopy to identify an unknown compound





Your notes



Your notes

Worked example

Look at the two infrared spectra below and determine which one corresponds to propanone and which one to propan-2-ol.



Copyright © Save My Exams. All Rights Reserved

Answer:

- IR spectrum **A** is **propanone** and spectrum **B** is **propan-2-ol**.
- In IR spectrum **A** the presence of a strong, sharp absorption around 1710 cm^{-1} corresponds to the characteristic $\text{C}=\text{O}$, carbonyl, group in a ketone.

- In spectrum **B** the presence of a strong, broad absorption around $3200\text{--}3500\text{ cm}^{-1}$ suggests that there is an alcohol group present, which corresponds to the -OH group in propan-2-ol.



Your notes